

HTATIP2 antibody

Product Information

Catalog No.:	FNab08714
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	HIV-1 Tat interactive protein 2, 30kDa
Synonyms:	Oxidoreductase HTATIP2 30 kDa HIV-1 TAT-interacting protein HIV-1 TAT-interactive protein 2 HTATIP2 CC3 TIP30
Observed MW:	30 kDa
Uniprot ID :	Q9BUP3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08714(TIP30 Antibody) at dilution of 1:600