

## TMCC1 antibody

### Product Information

Catalog No.:	FNab08735
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	transmembrane and coiled-coil domain family 1
Synonyms:	Transmembrane and coiled-coil domains protein 1 TMCC1 KIAA0779
Observed MW:	72 kDa
Uniprot ID :	O94876

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab08735(TMCC1 antibody) at dilution of 1:500