

## TMED3 antibody

### Product Information

Catalog No.:	FNab08739
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Potential role in vesicular protein trafficking, mainly in the early secretory pathway. Contributes to the coupled localization of TMED2 and TMED10 in the cis-Golgi network.

### Immunogen information

Immunogen:	transmembrane emp24 protein transport domain containing 3
Synonyms:	Transmembrane emp24 domain-containing protein 3 Membrane protein p24B p24 family protein gamma-4 (p24gamma4) p26 TMED3 C15orf22
Observed MW:	21-25 kDa
Uniprot ID :	Q9Y3Q3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab08739(TMED3 antibody) at dilution of 1:1000