

TMEM132A antibody

Product Information

Catalog No.:	FNab08751
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in embryonic and postnatal development of the brain. Increased resistance to cell death induced by serum starvation in cultured cells. Regulates cAMP-induced GFAP gene expression via STAT3 phosphorylation(By similarity).

Immunogen information

Immunogen:	transmembrane protein 132A
Synonyms:	Transmembrane protein 132A HSPA5-binding protein 1 TMEM132A HSPA5BP1 KIAA1583
Observed MW:	83-110 kDa
Uniprot ID :	Q24JP5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08751(TMEM132A Antibody) at dilution of 1:500