

TMEM165 antibody

Product Information

Catalog No.:	FNab08756
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as a calcium/proton transporter involved in calcium and in lysosomal pH homeostasis. Therefore, it may play an indirect role in protein glycosylation.

Immunogen information


Immunogen:	transmembrane protein 165
Synonyms:	Putative divalent cation/proton antiporter TMEM165 Transmembrane protein 165 Transmembrane protein PT27 Transmembrane protein TPARL TMEM165 TPARL
Observed MW:	35-40 kDa
Uniprot ID :	Q9HC07

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000; IF: 1:20-1:200
Image:	



IP Result of anti-TMEM165 (IP:FNab08756, 4ug;
Detection:FNab08756 1:1000) with HeLa cells
lysate 3440ug.



Jurkat cells were subjected to SDS PAGE followed
by western blot with FNab08756(TMEM165
antibody) at dilution of 1:3000