

TMEM176B antibody

Product Information

Catalog No.:	FNab08762
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in the process of maturation of dendritic cells. Required for the development of cerebellar granule cells(By similarity).

Immunogen information


Immunogen:	transmembrane protein 176B
Synonyms:	Transmembrane protein 176B Protein LR8 TMEM176B LR8
Observed MW:	30 kDa
Uniprot ID :	Q3YBM2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:200-1:1000; IP: 1:500-1:1000
Image:	



IP Result of anti-TMEM176B (IP: FNab08762, 4ug; Detection: FNab08762 1:500) with mouse lung tissue lysate 4000ug.



A549 cells were subjected to SDS PAGE followed by western blot with FNab08762(TMEM176B Antibody) at dilution of 1:500