

TMEM181 antibody

Product Information

Catalog No.:	FNab08763
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Mediates action of cytolethal distending toxins(CDT), which are secreted by many pathogenic bacteria. Expression level of TMEM181 is rate-limiting for intoxication.

Immunogen information


Immunogen:	transmembrane protein 181
Synonyms:	Transmembrane protein 181 TMEM181 GPR178 KIAA1423
Observed MW:	75 kDa
Uniprot ID :	Q9P2C4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab08763(TMEMP181 Antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab08763(TMEMP181 antibody) at dilution of 1:300