

AGMO antibody

Product Information

Catalog No.:	FNab08767
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	transmembrane protein 195
Synonyms:	Alkylglycerol monooxygenase Transmembrane protein 195 AGMO TMEM195
Observed MW:	52 kDa
Uniprot ID :	Q6ZNB7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:100-1:400
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab08767(TMEM195 antibody) at dilution of 1:200



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab08767(TMEM195 antibody) at dilution of 1:600