

## TMEM199 antibody

### Product Information

Catalog No.:	FNab08768
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in Golgi homeostasis.

### Immunogen information


Immunogen:	transmembrane protein 199
Synonyms:	Transmembrane protein 199 TMEM199 C17orf32
Observed MW:	18-23 kDa
Uniprot ID :	Q8N511

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab08768 (TMEM199 Antibody) at dilution of 1:50



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab08768 (TMEM199 Antibody) at dilution of 1:600