

TMEM208 antibody

Product Information

Catalog No.:	FNab08770
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as a negative regulator of endoplasmic reticulum-stress induced autophagy.

Immunogen information

Immunogen:	transmembrane protein 208
Synonyms:	Transmembrane protein 208 TMEM208
Observed MW:	19-29 kDa
Uniprot ID :	Q9BTX3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08770(TMEM208 Antibody) at dilution of 1:300