

TMEM214 antibody

Product Information

Catalog No.:	FNab08772
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Critical mediator, in cooperation with CASP4, of endoplasmic reticulum-stress induced apoptosis. Required for the activation of CASP4 following endoplasmic reticulum stress.

Immunogen information

Immunogen:	transmembrane protein 214
Synonyms:	Transmembrane protein 214 TMEM214
Observed MW:	77 kDa
Uniprot ID :	Q6NUQ4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



C2C12 cells were subjected to SDS PAGE followed by western blot with FNab08772(TMEM214 antibody) at dilution of 1:500