

TMEM50B antibody

Product Information

Catalog No.:	FNab08784
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	transmembrane protein 50B
Synonyms:	Transmembrane protein 50B HCV p7-trans-regulated protein 3 TMEM50B C21orf4
Observed MW:	14 kDa
Uniprot ID :	P56557

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08784(TMEM50B antibody) at dilution of 1:200