

## TMLHE antibody

### Product Information

Catalog No.:	FNab08798
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Converts trimethyllysine(TML) into hydroxytrimethyllysine(HTML)(PubMed:11431483, PubMed:23092983).

### Immunogen information

Immunogen:	trimethyllysine hydroxylase, epsilon
Synonyms:	Trimethyllysine dioxygenase, mitochondrial Epsilon-trimethyllysine 2-oxoglutarate dioxygenase Epsilon-trimethyllysine hydroxylase TML hydroxylase TML-alpha-ketoglutarate dioxygenase (TML dioxygenase, TMLD) TMLHE TMLH
Observed MW:	44 kDa
Uniprot ID :	Q9NVH6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human liver tissue were subjected to SDS PAGE followed by western blot with FNab08798(TMLHE antibody) at dilution of 1:500