

## anti- TMSB4X antibody

### Product Information

Catalog No.:	FNab08810
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Plays an important role in the organization of the cytoskeleton(By similarity). Binds to and sequesters actin monomers(G actin) and therefore inhibits actin polymerization. Seraspenide inhibits the entry of hematopoietic pluripotent stem cells into the S-phase.

### Immunogen information

Immunogen:	thymosin beta 4, X-linked
Synonyms:	TB4X, THYB4, TMSB4
Observed MW:	5 kDa
UniprotID :	P62328

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab08810(TMSB4X antibody) at dilution of 1:2000