

## TMTC1 antibody

### Product Information

Catalog No.:	FNab08811
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	transmembrane and tetratricopeptide repeat containing 1
Synonyms:	Protein O-mannosyl-transferase TMTC1 Transmembrane and TPR repeat-containing protein 1 TMTC1
Observed MW:	99 kDa, 65 kDa
Uniprot ID :	Q8IUR5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse bladder tissue were subjected to SDS PAGE followed by western blot with FNab08811(TMTC1 antibody) at dilution of 1:1000