

TNFRSF13B antibody

Product Information

Catalog No.:	FNab08825
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Receptor for TNFSF13/APRIL and TNFSF13B/TALL1/BAFF/BLYS that binds both ligands with similar high affinity. Mediates calcineurin-dependent activation of NF-AT, as well as activation of NF-kappa-B and AP-1. Involved in the stimulation of B-and T-cell function and the regulation of humoral immunity.

Immunogen information

Immunogen:	tumor necrosis factor receptor superfamily, member 13B
Synonyms:	Tumor necrosis factor receptor superfamily member 13B Transmembrane activator and CAML interactor TNFRSF13B TACI
Observed MW:	30-32 kDa
Uniprot ID :	O14836

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab08825(TNFRSF13B antibody) at dilution of 1:500