

TNP1 antibody

Product Information

Catalog No.:	FNab08843
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


In the elongating spermatids of mammals, the conversion of nucleosomal chromatin to the compact, non-nucleosomal form found in the sperm nucleus is associated with the appearance of a small set of basic chromosomal transition proteins.

Immunogen information

Immunogen:	transition protein 1(during histone to protamine replacement)
Synonyms:	Spermatid nuclear transition protein 1 (STP-1, TP-1) TNP1
Observed MW:	6-10 kDa
Uniprot ID :	P09430

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08843(TNP1 antibody) at dilution of 1:1000