

CELF3 antibody

Product Information

Catalog No.:	FNab08845
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	trinucleotide repeat containing 4
Synonyms:	CUGBP Elav-like family member 3 (CELF-3) Bruno-like protein 1 CAG repeat protein 4 CUG-BP- and ETR-3-like factor 3 ELAV-type RNA-binding protein 1 (ETR-1) Expanded repeat domain protein CAG/CTG 4 RNA-binding protein BRUNOL-1 Trinucleotide repeat-containing gene 4 protein CELF3 BRUNOL1 CAGH4 ERDA4 TNRC4
Observed MW:	55 kDa
Uniprot ID :	Q5SZQ8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cerebellum tissue slide using FNab08845(TNRC4 Antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab08845(TNRC4 Antibody) at dilution of 1:1000