

TOB1 antibody

Product Information

Catalog No.:	FNab08850
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Anti-proliferative protein; the function is mediated by association with deadenylase subunits of the CCR4-NOT complex.

Immunogen information


Immunogen:	transducer of ERBB2, 1
Synonyms:	Protein Tob1 Transducer of erbB-2 1 TOB1 TOB TROB1
Observed MW:	38-45 kDa
Uniprot ID :	P50616

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach using FNab08850(TOB1 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08850(TOB1 antibody) at dilution of 1:1000