

FNBP1L antibody

Product Information

Catalog No.:	FNab08852
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	formin binding protein 1-like
Synonyms:	Formin-binding protein 1-like Transducer of Cdc42-dependent actin assembly protein 1 (Toca-1) FNBP1L C1orf39 TOCA1
Observed MW:	63-75 kDa
Uniprot ID :	Q5T0N5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08852(Toca 1 antibody) at dilution of 1:300