

anti- TOP2A-Specific antibody

Product Information

Catalog No.:	FNab08868
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	topoisomerase(DNA) II alpha 170kDa
Synonyms:	None
Observed MW:	175 kDa
Uniprot ID :	P11388

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:3000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human spleen using FNab08868(TOP2A-Specific antibody) at dilution of 1:100



IP Result of anti-TOP2A-Specific (IP:FNab08868, 4ug; Detection:FNab08868 1:500) with HeLa cells lysate 2000ug.

HeLa cells were subjected to SDS PAGE followed by western blot with FNab08868(TOP2A-Specific antibody) at dilution of 1:400

