

TOR2A antibody

Product Information

Catalog No.:	FNab08873
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Salusins-alpha and-beta may be endocrine and/or paracrine factors able to increase intracellular calcium concentrations and induce cell mitogenesis. Salusins may also be potent hypotensive peptides.

Immunogen information

Immunogen:	torsin family 2, member A
Synonyms:	Prosalusin Torsin family 2 member A Torsin-2A Salusin-alpha Salusin-beta TOR2A
Observed MW:	40-45 kDa
Uniprot ID :	Q8N2E6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08873(TOR2A antibody) at dilution of 1:100