

TOR1A antibody

Product Information

Catalog No.:	FNab08876
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	torsin family 1, member A(torsin A)
Synonyms:	Torsin-1A Dystonia 1 protein Torsin ATPase-1A Torsin family 1 member A TOR1A DQ2 DYT1 TA TORA
Observed MW:	42-50 kDa
Uniprot ID :	O14656

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab08876(TOR1A antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08876(TOR1A Antibody) at dilution of 1:1000