

## TOX2 antibody

### Product Information

Catalog No.:	FNab08877
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Putative transcriptional activator involved in the hypothalamo-pituitary-gonadal system.

### Immunogen information


Immunogen:	TOX high mobility group box family member 2
Synonyms:	TOX high mobility group box family member 2 Granulosa cell HMG box protein 1 (GCX-1) TOX2 C20orf100 GCX1
Observed MW:	70 kDa
Uniprot ID :	Q96NM4

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-TOX2 (IP:FNab08877, 4ug;  
Detection:FNab08877 1:300) with mouse liver  
tissue lysate 4000ug.



mouse liver tissue were subjected to SDS PAGE  
followed by western blot with FNab08877(TOX2  
Antibody) at dilution of 1:300