

TPI1 antibody

Product Information

Catalog No.:	FNab08887
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2C10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	triosephosphate isomerase 1
Synonyms:	Triosephosphate isomerase (TIM) Methylglyoxal synthase Triose-phosphate isomerase TPI1 TPI
Observed MW:	27 kDa
Uniprot ID :	P60174

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:1000-1:4000
Image:	

HepG2 cells were subjected to SDS PAGE followed by western blot with FNab08887(TPI1 antibody) at dilution of 1:2000

