

TPPP antibody

Product Information

Catalog No.:	FNab08898
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in the polymerization of tubulin into microtubules, microtubule bundling and the stabilization of existing microtubules, thus maintaining the integrity of the microtubule network. May play a role in mitotic spindle assembly and nuclear envelope breakdown.

Immunogen information


Immunogen:	tubulin polymerization promoting protein
Synonyms:	Tubulin polymerization-promoting protein (TPPP) 25 kDa brain-specific protein TPPP/p25 p24 p25-alpha TPPP TPPP1
Observed MW:	24-28 kDa
Uniprot ID :	O94811

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab08898(TPPP Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab08898(TPPP Antibody) at dilution of 1:1000