

TPT1 antibody

Product Information

Catalog No.:	FNab08907
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in calcium binding and microtubule stabilization.

Immunogen information


Immunogen:	tumor protein, translationally-controlled 1
Synonyms:	Translationally-controlled tumor protein (TCTP) Fortilin Histamine-releasing factor (HRF) p23 TPT1
Observed MW:	22 kDa
Uniprot ID :	P13693

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab08907(TPT1 Antibody) at dilution of 1:200



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab08907(TPT1 antibody) at dilution of 1:1000