

TRA2A antibody

Product Information

Catalog No.:	FNab08913
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Sequence-specific RNA-binding protein which participates in the control of pre-mRNA splicing.

Immunogen information

Immunogen:	transformer 2 alpha homolog(Drosophila)
Synonyms:	Transformer-2 protein homolog alpha (TRA-2 alpha, TRA2-alpha) Transformer-2 protein homolog A TRA2A
Observed MW:	33 kDa
Uniprot ID :	Q13595

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab08913(TRA2A antibody) at dilution of 1:3000