

## TRAPPC4 antibody

### Product Information

Catalog No.:	FNab08948
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in vesicular transport from endoplasmic reticulum to Golgi.

### Immunogen information

Immunogen:	trafficking protein particle complex 4
Synonyms:	Trafficking protein particle complex subunit 4 Hematopoietic stem/progenitor cell protein 172 Synbindin TRS23 homolog TRAPPC4 SBDN
Observed MW:	20 kDa
Uniprot ID :	Q9Y296

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08948(TRAPPC4 antibody) at dilution of 1:400