

TRAPPC9 antibody

Product Information

Catalog No.:	FNab08950
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	trafficking protein particle complex 9
Synonyms:	Trafficking protein particle complex subunit 9 NIK- and IKBKB-binding protein Tularik gene 1 protein TRAPPC9 KIAA1882 NIBP
Observed MW:	128 kDa, 139 kDa
Uniprot ID :	Q96Q05

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Neuro-2a cells were subjected to SDS PAGE followed by western blot with FNab08950(Trappc9,NIBP antibody) at dilution of 1:300