

TRIB2 antibody

Product Information

Catalog No.:	FNab08965
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Interacts with MAPK kinases and regulates activation of MAP kinases. Does not display kinase activity(By similarity).

Immunogen information

Immunogen:	tribbles homolog 2(Drosophila)
Synonyms:	Tribbles homolog 2 (TRB-2) TRIB2 TRB2
Observed MW:	42 kDa
Uniprot ID :	Q92519

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08965 (TRIB2 antibody) at dilution of 1:500