

TRIM17 antibody

Product Information

Catalog No.:	FNab08970
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as a ubiquitin E3 ligase.

Immunogen information

Immunogen:	tripartite motif-containing 17
Synonyms:	E3 ubiquitin-protein ligase TRIM17 RING finger protein 16 RING-type E3 ubiquitin transferase TRIM17 Testis RING finger protein Tripartite motif-containing protein 17 TRIM17 RBCC RNF16 TERF
Observed MW:	52 kDa
Uniprot ID :	Q9Y577

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08970(TRIM17 antibody) at dilution of 1:300