

TRIM2 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab08971 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


UBE2D1-dependent E3 ubiquitin-protein ligase that mediates the ubiquitination of NEFL and of phosphorylated BCL2L11. Plays a neuroprotective function. May play a role in neuronal rapid ischemic tolerance.

Immunogen information


| | |
|--------------|--|
| Immunogen: | tripartite motif-containing 2 |
| Synonyms: | Tripartite motif-containing protein 2 E3 ubiquitin-protein ligase TRIM2 RING finger protein 86 RING-type E3 ubiquitin transferase TRIM2 TRIM2 KIAA0517 RNF86 |
| Observed MW: | 70-82 kDa |
| Uniprot ID : | Q9C040 |

Application

| | |
|-----------------------|---|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IF, IHC |
| Recommended dilution: | WB: 1:1000-1:10000; IP: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100 |
| Image: | |



Immunohistochemistry of paraffin-embedded human kidney slide using FNab08971 (TRIM2 Antibody) at dilution of 1:50



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with FNab08971 (TRIM2 antibody) at dilution of 1:500