

## TRIM31 antibody

### Product Information

Catalog No.:	FNab08979
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Regulator of Src-induced anchorage independent cell growth(By similarity). May have E3 ubiquitin-protein ligase activity.

### Immunogen information


Immunogen:	tripartite motif-containing 31
Synonyms:	E3 ubiquitin-protein ligase TRIM31 Tripartite motif-containing protein 31 TRIM31 C6orf13 RNF
Observed MW:	48 kDa
Uniprot ID :	Q9BZY9

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab08979(TRIM31 antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab08979(TRIM31 antibody) at dilution of 1:400