

TRIM38 antibody

Product Information

Catalog No.:	FNab08983
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


E3 ubiquitin-protein ligase. Mediates 'Lys-48'-linked polyubiquitination and proteasomal degradation of the critical TLR adapter TICAM1, inhibiting TLR3-mediated type I interferon signaling.

Immunogen information

Immunogen:	tripartite motif-containing 38
Synonyms:	E3 ubiquitin-protein ligase TRIM38 RING finger protein 15 Tripartite motif-containing protein 38 Zinc finger protein RoRet TRIM38 RNF15 RORET
Observed MW:	53 kDa
Uniprot ID :	O00635

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000; IF: 1:50-1:500
Image:	



human liver tissue were subjected to SDS PAGE followed by western blot with FNab08983(TRIM38 antibody) at dilution of 1:800