

anti- TRIM69 antibody

Product Information

Catalog No.:	FNab08994
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May have E3 ubiquitin-protein ligase activity. May play a role in apoptosis.

Immunogen information

Immunogen:	tripartite motif-containing 69
Synonyms:	HSD 34, HSD34, RING finger protein 36, RNF36, Trif, TRIM69, tripartite motif containing 69
Observed MW:	60-70 kDa
UniprotID :	Q86WT6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human gliomas using FNab08994(TRIM69 antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08994(TRIM69 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com