

NTRK3 antibody

Product Information

Catalog No.:	FNab09003
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	neurotrophic tyrosine kinase, receptor, type 3
Synonyms:	NT-3 growth factor receptor GP145-TrkC (Trk-C) Neurotrophic tyrosine kinase receptor type 3 TrkC tyrosine kinase NTRK3 TRKC
Observed MW:	92 kDa
Uniprot ID :	Q16288

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



fetal human brain tissue were subjected to SDS
PAGE followed by western blot with
FNab09003(NTRK3 Antibody) at dilution of 1:300