

## TRNAU1AP antibody

### Product Information

Catalog No.:	FNab09010
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	tRNA selenocysteine 1 associated protein 1
Synonyms:	tRNA selenocysteine 1-associated protein 1 SECP43 tRNA selenocysteine-associated protein 1 TRNAU1AP SECP43 TRSPAP1
Observed MW:	36 kDa
Uniprot ID :	Q9NX07

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab09010(TRNAU1AP,SECP43 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09010(TRNAU1AP,SECP43 antibody) at dilution of 1:1500