

TROAP antibody

Product Information

Catalog No.:	FNab09011
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Could be involved with bystin and trophinin in a cell adhesion molecule complex that mediates an initial attachment of the blastocyst to uterine epithelial cells at the time of the embryo implantation.

Immunogen information

Immunogen:	trophinin associated protein(tastin)
Synonyms:	Tastin Trophinin-assisting protein Trophinin-associated protein TROAP
Observed MW:	110 kDa
Uniprot ID :	Q12815

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:10-1:100
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09011(TROAP Antibody) at dilution of 1:300