

TRPT1 antibody

Product Information

Catalog No.:	FNab09027
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the last step of tRNA splicing, the transfer of the splice junction 2'-phosphate from ligated tRNA to NAD to produce ADP-ribose 1"-2" cyclic phosphate.

Immunogen information

Immunogen:	tRNA phosphotransferase 1
Synonyms:	tRNA 2'-phosphotransferase 1 TRPT1
Observed MW:	None
Uniprot ID :	Q86TN4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab09027 (TRPT1 Antibody) at dilution of 1:50