

## TRUB2 antibody

### Product Information

Catalog No.:	FNab09032
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be responsible for synthesis of pseudouridine from uracil in transfer RNAs.

### Immunogen information


Immunogen:	TruB pseudouridine(psi) synthase homolog 2(E. coli)
Synonyms:	Pseudouridylate synthase TRUB2, mitochondrial TruB pseudouridine synthase homolog 2 tRNA pseudouridine 55 synthase TRUB2 (Psi55 synthase TRUB2) TRUB2
Observed MW:	37 kDa
Uniprot ID :	O95900

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin using FNab09032(TRUB2 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09032(TRUB2 antibody) at dilution of 1:500