

TSFM antibody

Product Information

Catalog No.:	FNab09039
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Associates with the EF-Tu.GDP complex and induces the exchange of GDP to GTP. It remains bound to the aminoacyl-tRNA.EF-Tu.GTP complex up to the GTP hydrolysis stage on the ribosome.

Immunogen information


Immunogen:	Ts translation elongation factor, mitochondrial
Synonyms:	Elongation factor Ts, mitochondrial (EF-Ts, EF-TsMt) TSFM
Observed MW:	39 kDa
Uniprot ID :	P43897

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNa09039(TSFM antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNa09039(TSFM Antibody) at dilution of 1:600