

TSNAX antibody

Product Information

Catalog No.:	FNab09049
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts in combination with TSN as an endonuclease involved in the activation of the RNA-induced silencing complex(RISC). Possible role in spermatogenesis.

Immunogen information

Immunogen:	translin-associated factor X
Synonyms:	Translin-associated protein X Translin-associated factor X TSNAX TRAX
Observed MW:	77 kDa
Uniprot ID :	Q99598

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200;IF: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human gliomas using FNab09049(TSNAX antibody) at dilution of 1:50