

TSPAN33 antibody

Product Information

Catalog No.:	FNab09063
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays an important role in normal erythropoiesis. It has a role in the differentiation of erythroid progenitors(By similarity).

Immunogen information


Immunogen:	tetraspanin 33
Synonyms:	Tetraspanin-33 (Tspan-33) Penumbra (hPen) Proerythroblast new membrane TSPAN33 PEN
Observed MW:	70 kDa
Uniprot ID :	Q86UF1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab09063(TSPAN33 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab09063(TSPAN33 antibody) at dilution of 1:400