

TSPAN6 antibody

Product Information

Catalog No.:	FNab09065
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	tetraspanin 6
Synonyms:	Tetraspanin-6 (Tspan-6) A15 homolog Putative NF-kappa-B-activating protein 321 T245 protein Tetraspanin TM4-D Transmembrane 4 superfamily member 6 TSPAN6 TM4SF6
Observed MW:	28 kDa
Uniprot ID :	O43657

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab09065(TSPAN6 antibody) at dilution of 1:1000