

TSPYL1 antibody

Product Information

Catalog No.:	FNab09068
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	TSPY-like 1
Synonyms:	Testis-specific Y-encoded-like protein 1 (TSPY-like protein 1) TSPYL1 TSPYL
Observed MW:	49 kDa
Uniprot ID :	Q9H0U9

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab09068(TSPYL1 antibody) at dilution of 1:500