

## TTC31 antibody

### Product Information

Catalog No.:	FNab09081
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	tetratricopeptide repeat domain 31
Synonyms:	Tetratricopeptide repeat protein 31 (TPR repeat protein 31) TTC31
Observed MW:	57 kDa
Uniprot ID :	Q49AM3

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab09081(TTC31 antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab09081(TTC31 antibody) at dilution of 1:1000