

TTC39A antibody

Product Information

Catalog No.:	FNab09085
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	tetratricopeptide repeat domain 39A
Synonyms:	Tetratricopeptide repeat protein 39A (TPR repeat protein 39A) Differentially expressed in MCF-7 with estradiol protein 6 (DEME-6) TTC39A C1orf34 KIAA0452
Observed MW:	66 kDa
Uniprot ID :	Q5SRH9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab09085(TTC39A antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab09085(TTC39A antibody) at dilution of 1:100