

TTF1 antibody

Product Information

Catalog No.:	FNab09090
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4D3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Multifunctional nucleolar protein that terminates ribosomal gene transcription, mediates replication fork arrest and regulates RNA polymerase I transcription on chromatin. Plays a dual role in rDNA regulation, being involved in both activation and silencing of rDNA transcription. Interaction with BAZ2A/TIP5 recovers DNA-binding activity.

Immunogen information


Immunogen:	NK2 homeobox 1
Synonyms:	Transcription termination factor 1 (TTF-1) RNA polymerase I termination factor Transcription termination factor I (TTF-I) TTF1
Observed MW:	42 kDa
Uniprot ID :	P43699

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab09090(NKX2-1 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09090(NKX2-1 antibody) at dilution of 1:1000