

TTLL10 antibody

Product Information

Catalog No.:	FNab09097
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inactive polyglycyclase.

Immunogen information

Immunogen:	tubulin tyrosine ligase-like family, member 10
Synonyms:	Inactive polyglycyclase TTLL10 Tubulin--tyrosine ligase-like protein 10 TTLL10
Observed MW:	49 kDa
Uniprot ID :	Q6ZVT0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab09097(TTLL10 Antibody) at dilution of 1:1000